

# Sling Guide



## **BEAUCARE SLINGS**

### **OUR SLING RANGE**

Our sling range focuses on comfort, durability and fit. This guide is designed to aid you to choose the most suitable sling for your needs and environment.

- Choose between 10 designs,
   3 fabric choices and 6 sizes
- Colour-coded edges and loops
- Both loop and clip fixings available
- Embroidered labels that last the life of the sling



Our team
can help you
select the right
sling for your
needs

# CUSTOMERS LOVE OUR



Nationwide Fleet Delivery



Diverse Product Range



Dedicated Account Managers



Award Winning Service



Vicki Crombleholme
Training and Compliance Manager

# **ATTACHMENTS & FABRICS**

### **ATTACHMENT TYPES**



#### **LOOPS**

Various coloured loop straps provide adjustment and correct positioning angle for the user.



#### **CLIPS**

Should your hoist attach using clip fixings we are able to have compatible slings made to order.



### **FABRICS**



- Durable, low friction fabric
- Ideal for general transfers
- Withstands high laundry temperatures



- For use when bathing & showering
- Quickly drains water, facilitating fast drying



- Designed to be used and left in place between transfers
- Withstands high laundry temperatures

## **SLING RANGE**



#### UNIVERSAL SLING

A standard quick-fit, universal sling that works for 85% of individuals. Offers a more general solution to individuals lifting needs.



#### UNIVERSAL DELUXE SLING

Featuring a head support and padding on the leg area, this sling offers a higher level of comfort and support than a standard universal sling.



#### CLASSIC HAMMOCK SLING

Designed to provide support and comfort by reducing hip movement.



#### **DELUXE HAMMOCK SLING**

A hybrid sling designed to offer the full support of a hammock sling with the function of a universal. It can be applied and removed while seated.




#### HAMMOCK SLING

Hammock style long seat sling. Designed for amputees and individuals who are not comfortable in a divided leg sling.

#### **OPTIONS AVAILABLE:**















#### IN CHAIR HAMMOCK SLING

Provides a secure fit with safety pommel. Suitable for those with amputations.

#### **OPTIONS AVAILABLE:**





#### IN CHAIR SPLIT LEG HAMMOCK SLING


Split leg design allowing you to fold away the leg supports preventing build up of heat and adding to user comfort.





#### **TOILETING SLING**


Designed with a larger aperture which allows for easy access during toileting. Clipped belt with a fastening at the waist.





#### TRANSPORTING SLING


Designed to offer greater support for an individual using a stand aid hoist. Clipped belt with a fastening at the waist.





#### **STAND AID SLING**

Designed to be used in conjunction with a stand aid hoist. Clipped belt with a fastening at the waist.



# **SIZE GUIDE**




### **USE THE SIZE TABLE BELOW FOR:**

Toileting Sling Transporting Sling Stand Aid Sling

SLING SIZE	USER HEIGHT*	USER WIDTH*
X SMALL	145-155	70-85
SMALL	150-170	75-90
MEDIUM	160-180	85-100
LARGE	160-185	90-115
X LARGE	170-195	100-125






#### **USE THE SIZE TABLE BELOW FOR:**

Universal Sling
Deluxe Universal Sling
Classic Hammock Sling
Deluxe Hammock Sling
In Chair Hammock Sling
In Chair Split Leg Hammock Sling

SLING SIZE	SLING WIDTH*	USER BACK*	USER WIDTH*
X SMALL	77	49-52	45-54
SMALL	86	56-58	54-60
MEDIUM	96	62-68	60-72
LARGE	116	67-73	72-82
X LARGE	129	70-73	82-100

<sup>\*</sup>All measurements are in cm




#### **USE THE SIZE TABLE BELOW FOR:**

Hammock Sling

SLING SIZE	SLING WIDTH*	USER BACK*	USER WIDTH*
X SMALL	81	104	45-54
SMALL	86	109	54-60
MEDIUM	96	114	60-68
LARGE	106	122	68-76
X LARGE	116	129	76-92

<sup>\*</sup>All measurements are in cm







Beaucare® Medical Ltd Hornbeam Square West, Hornbeam Park, Harrogate, North Yorkshire, HG2 8PA

> t: 01423 873666 e: sales@beaucare.com w: www.beaucare.com